

Quick Reference Guide

1. Coding Your Meter

The display will show "no CODE" the first time batteries are installed.

⚠ The meter is for outside the body (in vitro diagnostic use only) use.

Insert the Code Card into the slot. The screen displays "OK".

⚠ Please read the User's Manual and Test Strip instructions before test.

The code number will be shown on the display.

Make sure the code number on the display matches the code number on the test strip bottle.

2. Running A Control Solution Test

Make sure the code number on the Code Card matches the code number on the test strip bottle.

Insert the test strip into the meter.

⚠ Replace the bottle cap immediately after removing a test strip from the bottle.

Squeeze a tiny drop on the strip's reaction zone.

Check if the reading falls within the range printed on the test strip bottle.

3. Inserting A Lancet into the Lancing Device

Take the cap off.

⚠ To disinfect the cap after cleaning by placing it in 70% rubbing alcohol for 10 minutes. Allow the cap to air-dry after disinfecting.

Insert the lancet into the lancing device firmly.

Twist off the protective cover.



Please read the User's Manual and Test Strip Instructions before testing.

If you have any questions and/or need assistance, please contact your healthcare professional.

Select the desired penetration depth.

Pull on the sliding barrel of the lancing device until it clicks and then release.

4. Alternative Site Testing

This graphic shows where the meter cleared for alternative site testing.

⚠ Check with your healthcare professional before test. other sites that not fingertip.

Massage the puncture area of forearm or palm for a few minutes.

Press and hold the lancing device with one clear adjustable tip against the forearm.

If you choose the palm, repeat the preceding three steps to obtain a drop of blood sample.

Press the trigger button to activate the lancing device.

⚠ Please do NOT use the first drop of blood sample.

Hold the device against forearm or palm and increase pressure until the blood sample size is sufficient.

5. Testing Your Blood Glucose

Insert the test strip into the meter.

⚠ For accurate results, wash your hands with soap and warm water. Rinse and dry thoroughly.

Press the trigger button to activate the lancing device.

⚠ Replace the bottle cap immediately after removing a test strip from the bottle.

Be sure to get enough blood on the strip's reaction zone.

The result will appear on the display after 5 seconds.

